



Value Life

Introducing...
smartmidline™

CT-rated midline
Up to 325psi and 7ml/s*

New



Introducing smartmidline™

CT-rated midline Up to 325psi and 7ml/s*

Easy and quick to place

Direct Seldinger technique

Midline label

For easy identification

Suitable for all patient and user needs

From 2 to 5Fr and 4 to 25cm

Atraumatic insertion


Rounded or conical flexible tip, reduces the risk of thrombosis

Integrated extension line

Keeps manipulation away from the insertion site, prevents infection risks.



*Only 5Fr catheters



Kit composition*

- Puncture needle, 21G - 4.5cm
- Stainless steel straight guidewire, 0.018" (0.46mm)
- Dilator (only included with 4Fr catheter)

*5Fr contains catheter only.

Implement smartmidline today

Contact us today to see how smartmidline can help you

Call: 01793 748830 or **Email:** marketing@vygon.co.uk

Product Information

Vygon Code	Catheter Size	Length (cm)	Guidewire Length (cm)	Max. Flowrate (Viscosity 11.8 cP)
128.1204	2Fr	4	23	1.5ml/s
128.1206		6	23	1.5ml/s
128.12062		6	40	1.5ml/s
128.1208		8	26	1ml/s
128.12082		8	40	1ml/s
128.1210		10	26	1ml/s
128.1215		15	40	0.5ml/s
128.1220		20	50	0.5ml/s
Coming soon	3Fr	6	23	TBC
Coming soon		8	26	
Coming soon		10	26	
Coming soon		15	40	
Coming soon		20	TBC	
Coming soon		25	TBC	
Coming soon	4Fr	12	40	TBC
128.1415		15	50	5ml/s
128.1420		20	50	5ml/s
Coming soon		25	60	TBC
128.1515	5Fr*	15	-	7ml/s
128.1520		20	-	7ml/s

*5Fr catheters are delivered without accessories. Must be used with an MST kit with a 5.5Fr introducer sheath.

For further information, please contact: marketing@vygon.co.uk

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon (UK) Ltd, The Pierre Simonet Building, V Park, Gateway North,
Latham Road, Swindon, Wiltshire SN25 4DL

Tel: 01793 748800 Fax: 01793 748899 Email: vygon@vygon.co.uk

vygon.co.uk [@vygonuk](https://twitter.com/vygonuk) [in vygonuk](https://www.linkedin.com/company/vygonuk) [vygonuk](https://www.youtube.com/channel/UCvYgonuk) [vygonuk](https://www.instagram.com/vygonuk)

Copyright Vygon (UK) Ltd 2017 Content correct as of: 12/2017 Code: PS463 v1



Value Life