

Safe disconnection of an Elastomeric Pump in Homecare

Context

Elastomeric pumps are frequently used for ambulatory chemotherapy treatments. This device presents a lot of advantages: time saving for the technician (pre-programmed device), mobility and comfort for the patient. But handling cytotoxic drugs represents a risk of chemical contamination, especially at patient's home, where no hospital protocols can be applied.



Problem

High risk of chemical contamination on disconnection for the nurse, the patient, their relatives and the environment.



Solution

Qimono connection prevents any contamination and preserves a safe environment for homecare treatments.



qimono♂



qimono♀



Safe disconnection

Ordering code

Vygon code	Description	Indication	Qty
7210.04	Set with Qimoprime, Qimo male ♂ and Qimo female ♀ in the same packaging	Can be used to secure the injection site of the elastomeric pump and the distal end of its tubing, but other applications are possible	50 box 600 case